



## Specimen Requirements

### Endometriosis Panel

TESTS	SPECIMEN	VACUTAINER TUBE	STORAGE TEMP	SHIPPING TEMP
<ul style="list-style-type: none"> <li>• Anti-phospholipid Antibodies (APA)</li> <li>• Anti-nuclear Antibodies screen (ANA)</li> <li>• Cancer Antigen 125 (CA 125)</li> </ul>	Serum	<b>Tiger SST</b> or <b>Red</b> spin & separate	Refrigerate 2-6° C	Ambient