

Specimen Requirements

FDA Mandated Donor Screening Panels

TESTS	SPECIMEN	VACUTAINER TUBE	STORAGE TEMP	SHIPPING TEMP
<ul style="list-style-type: none"> • HIV I Nucleic Acid Testing (NAT) • HCV (NAT) • HBV (NAT) • <i>WNV NAT</i> 	Whole Blood	Lavender Top	Refrigerate 2-6° C	Ambient
<ul style="list-style-type: none"> • Chlamydia Trachomatis NAT • Neisseria Gonorrhoea NAT 	Vaginal Swab Urine	Aptima Swab Pack Note: leave swab in transport media	Refrigerate 2-6° C	Refrigerate 2-6° C
<ul style="list-style-type: none"> • HIV1/2 plus O antibodies (anti-HIV 1/2 plus O) • Hepatitis B surface antigen (anti-HBcAg) • Hepatitis B Core Antigen (anti-HBcAg) • Hepatitis C Antibody (anti-HCV) • Trep Sure • HTLV I/II (anti-HTLV I/II) ♂ • Anti Cytomegalovirus (CMV) Total ♂ 	Serum	Tiger SST or Red spin & separate	Refrigerate 2-6° C	Ambient